

Stylix

Radiographic Table

The **STYLIX** radiographic table provides exceptional versatility for a wide range of x-ray procedures for busy radiology departments.

The compact space-saving design allows installation in even small x-ray rooms. **STYLIX**, with its large tabletop and extended travel, enables full radiographic coverage with minimal patient movement.

A two-step foot treadle provides maximum operator flexibility in controlling table movement. The four-way floating tabletop glides on a precision bearing system with electromagnetic locks to secure the tabletop and patient for quick and accurate positioning.

The flat tabletop permits easy patient transfer and unobstructed cross-table examinations. Durable and easy to use, this table is suitable for all environments from hospitals to private practices.

The **STYLIX** table is compatible with most radiographic systems and is designed to handle the rigors of the heaviest caseloads.

Control-X Medical offers revolutionary design, quality and value.



Imaging for the Future

The **STYLIX** radiographic table provides full radiographic coverage for maximum imaging flexibility.

Features:

Convenient control

A two-step foot treadle provides maximum flexibility for the technologist to achieve accurate patient positioning. Automatic center stop at radiographic position.

Tabletop movement

The large easy-to-clean tabletop supports up to 500 pounds while its precision bearing system provides effortless patient positioning. Combined with the extended Bucky travel, STYLIX allows for full coverage of the entire length of the table base.

Compatibility

The extended tabletop travel provides access for wheelchair patients and allows the use of C-Arms. It is compatible with most floor- and ceiling-mounted tube supports and so with most radiographic systems.

Construction

The low tabletop absorption material reduces patient exposure while the small skin-film distance permits precise low magnification film studies. Standard side rails support a wide variety of accessories.

Model diversity

Accepts size sensing Bucky's with automatic cassette centering feature, size sensing or non-size sensing reciprocating Bucky and tray system or a grid cabinet. Each with room for AEC chamber and grid.

Specifications:

Physical Dimensions:

Tabletop length:	86" (220 cm)
Optional short version:	73" (185 cm)
Tabletop width:	32" (81 cm)
Tabletop height:	29 1/4" (74.5 cm)
Film tabletop distance:	3" (7.5 cm)

Travel:

Tabletop:	
Longitudinal:	43" (110 cm)
Horizontal:	± 4 3/4" (12 cm)
Bucky:	22 1/2" (57 cm)

Cassette sizes:

All sizes up to a format of:	17"x17" (43 cm)
------------------------------	-----------------

Tabletop absorption density: Less than 0.85 mm Al equivalent

Power Requirements:

Voltage:	24 VAC ±10%
Frequency:	50/60 Hz
Current consumption:	1.5 A

Optional accessories:

Head clamp
Hip clamp
Compression belt
Lateral cassette holder



Specifications subject to change without notice

Your Dealer:

Control-X Medical Inc.

2289 Westbrooke Drive, Columbus, OH 43228

Phone: 614-777-9729 Fax: 614-777-0395

sales@cxmed.com Visit us on the Internet: www.cxmed.com