

VETELEV ELEVATING TABLE

A SIMPLIFIED AND COST-CONSCIOUS VERSION OF OUR POPULAR PHOENIX MEDICAL PATIENT TABLES, FEATURING



SUBSTANTIAL FREEDOM OF POSITIONING

- · 36 cm vertical movement
- Floating table top with electromagnetic brakes (152 cm x 71 cm with extension option)
- Over 100 kg lifting capacity





TOTAL SAFETY AND COMFORT

- Integrated crash guard sensors on the table • Grid cabinet for standard cassette size wifi or
- fixed DR panels ensuring crystal clear images · Configuration for fixed or removable grids, including safe "parking slot" for unused grid
- Wireless foot pedal (with long-life battery) to minimize cabling (optional)
- Foot pedal doubles as light field switch



ZOOMAX TUBE STAND

- · Sleek, robust design
- Fixed position tube with 100-136 cm SID (with Vetelev table)
- · Free-standing setup for total positioning freedom (optional)





INTEGRATED GENERATOR CONTROL AND IMAGE ACOUISITION

Single console seamless workflow:

- Unlimited APRs
- Full DICOM compatibility

Sophisticated image processing features:

- Detector status indicators
- · Wifi signal strength indicator
- · Wide variety of DR detector choices
- Portable solutions with tablet. laptop workstation
- Mini PACS included
- · Cloud storage available
- · Mobile viewing on any device

HIGH QUALITY COMPONENTS

- Industry-leading x-ray tubes
- High-density grids
- · Certified manual collimators with light beam or laser for easy positioning
- · Powerful generators for minimal exposure times

SELECTION OF CONTROL CONSOLES

- · 22" Touchscreen Monitor Console to set exposure parameters and for full image capture and analysis
- 10" LCD console integrated into X-ray tube cradle - to set exposure parameters, featuring:
 - slim, modern User Interface • customizable main page and workflow
 - buttons (e.g. APR, grid control, focus, etc.)
 - remote diagnostics and software updates easy installation
- · ZooMax LC or CPI Rad membrane consoles with comfort of all key generator parameter settings, including APRs and advanced controls such as AEC









GENERATOR OPTIONS

- · 32 kW to 50 kW output range
- Stored energy generators up to 50 kW (capacitor discharge) at 15 Amp
- 2.8 kW to 5 kW monoblock generators



Service software allows troubleshooting to minimize downtime.



Ready to connect (secure connection

CX Remote Supp...

Your ID

Password ++++



Please contact us at sales@cxmed.com or at +36 1 381 0301.

CERTIFICATIONS



Zoomax Vetelev Gold Info sheet - Version 2.0 November 2021

Patient positioning quide